# Veriti™

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
The Power of MORE

**AB** Applied Biosystems

**Veriti** 96 Well Fast Thermal Cycler

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
**Figure 1a.** Screenshot from Veriti™ Thermal Cycler showing the annealing temperature being set for 6 different primers.

**Figure 1b.** VeriFlex™ Blocks: Six individual peltier blocks.

**Figure 1c.** PCR results showing six primer sets run in a single PCR amplification cycle. Results indicate that the Veriti™ Thermal Cycler can run six assays at six annealing temperatures during the same PCR run.

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
Figure 3. Screenshot of default methods stored on the Veriti™ Thermal Cycler. Select a default method for quick setup. Default methods can also be edited to your requirements and saved for later use.



Figure 4. Edit run screen from Veriti™ Thermal Cycler. Simply touch a step and enter new values into keypad (not shown) to create new methods or to alter existing ones.

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
**Figure 5:** Gel image for p53 (1.4kb) comparing results from both fast and standard PCR amplification generated on the 0.2 mL Veriti™ 96-Well Thermal Cycler.

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**Figure 6:** Target p53 (1.4kb) generated using the Fast PCR Master Mix in 35 minutes versus 2 hours using standard PCR chemistry. All results were generated on the 0.2 mL Veriti 96-Well Thermal Cycler.

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
Description	Quantity	P/N
Applied Biosystems Veriti™ 96-Well Fast Thermal Cyclers, 0.1 mL	1 instrument	4375305
Applied Biosystems Veriti™ 96-Well Thermal Cyclers, 0.2 mL	1 instrument	4375786
Applied Biosystems Veriti™ 384-Well Thermal Cyclers, 0.02 mL	1 instrument	4388444
Applied Biosystems Veriti™ 60-Well, 0.5mL Thermal Cyclers	1 instrument	4384638
9 Channel Temperature Verification Kit, for VeriFlex™ 0.1 mL 96 Well Sample Block	1 kit	4373997
9 Channel Temperature Verification Probe, for VeriFlex™ 0.1 mL 96 Well Sample Block, Upgrade Kit	1 kit	4375938
9 Channel Temperature Verification Kit, for VeriFlex™ 0.2 mL 96 Well Sample Block	1 kit	4377669
9 Channel Temperature Verification Probe, for VeriFlex™ 0.2 mL 96 Well Sample Block, Upgrade Kit	1 kit	4377670
9 Channel Temperature Verification Kit for 384-Well, 0.02 mL Sample Block	1 kit	4385338
9 Channel Temperature Verification Probe for 384-Well, 0.02 mL Sample Block, Upgrade Kit	1 kit	4385339
60-Well, 0.5 mL Temperature Verification Kit	1 kit	4309924

如需了解更多Veriti™ 热循环仪的信息  
请访问 [info.appliedbiosystems.com/veriti](http://info.appliedbiosystems.com/veriti)

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