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Run Method: 9700-Max-Mode

Folder: Default

Reaction Volume: 10 Cover Temperature: 105.0

Notes:

Save & Exit

Enter a name, reaction volume, and choose a folder for this run method. Notes can also be added. When you are finished, touch 'Save & Exit'.

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Tools Menu

View Last Run Calculate T_m Convert a Method Run TNU Test

Run Temperature Verification Run Heated Cover Verification Run Cycle Performance Test Show Statistics

Calibrate Touch Screen To Previous Menu

Log Off 2007-02-14 | 8:22 PM | User: gregor

B. 点击“Calculate T_m”

T_m Calculator

Salt Concentration: 50 mM Primer Concentration: 0.20 uM

Primer 1 Sequence:

Primer 2 Sequence:

Calculate

Touch the fields to enter concentrations and the primer sequences, then touch Calculate T_m to calculate the melting temperatures. Maximum input for each primer sequence is 30.

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The screenshot shows the 'Tm Calculator' interface. At the top, there is a home icon and a close button (X). Below this, there are two input fields: 'Salt Concentration: 50 mM' and 'Primer Concentration: 0.20 uM'. Underneath are two text input fields for primer sequences: 'Primer 1 Sequence: CGTTTAGCT' and 'Primer 2 Sequence: CCCCCCGTATATTTTGC'. At the bottom of the input area, it displays 'TM for Primer 1: +17.6' and 'TM for Primer 2: +59.8', followed by a yellow 'Calculate' button. A green question mark icon is located at the bottom right. Below the input area, there is a small text block: 'Touch the fields to enter concentrations and the primer sequences, then touch Calculate Tm to calculate the melting temperatures. Maximum input for each primer sequence is 30.'