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
Select your protocol of interest from the preloaded list.

Versatility and Performance

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
Easily read multiplexed cell viability, apoptosis and cytotoxicity assays.

Multiplexed cytotoxicity and viability assays performed on digitonin-treated K562 cells.

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Efficiently read the output of kinase inhibitor screening and automatically analyze results.


Receptor tyrosine kinase inhibition was assessed following incubation with gefitinib, a selective inhibitor of EGFR, using the ADP-Glo™ Kinase Assay.

State-of-the-Art Sensitivity

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
Detect luminescence from Promega bioassays and evaluate mAb activity.

ADCC Bioassay Target cells were incubated with a series of concentrations of Control Ab, Anti-CD20, followed by addition of ADCC Bioassay Effector Cells.

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Analyze the interaction of proteins using BRET.

A dose-response curve showing the disruption of the p53 and MDM2 interaction using the inhibitor Nutlin-3 in a NanoBRET™ assay.

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SiLA Rapid Integration®

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Can be integrated with:

- Tecan Freedom EVO®
- Hudson Robotics Solo™
- PlateCraneEX™
- Hamilton Microlab® STAR series
- many other hardware systems

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Luminescence Module Specifications

Detector	Top-reading, head-on photon counting photomultiplier tube (PMT)
Wavelength Range	350–700nm
Detection Limit	3×10^{-21} moles of luciferase
Linear Dynamic Range	9 logs
Cross talk	Less than 3×10^{-5} (white, 96-well plate, Corning 3912)
Filtered Luminescence	GloMax [®] Discover: built-in filter paddle
Dynamic Range	9 logs
Sensitivity	3×10^{-21} moles of luciferase

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Fluorescence Module Specifications

Detector	PIN-photodiode
Light Source	Wavelength-matched LED
Read Position	Top reading
Wavelengths selection	Filter module with 5 standard excitation and emission filters. Empty filter positions are available for custom excitation and emission wavelengths.
Wavelengths included	UV (Ex: 365nm, Em: 415–445nm) Blue (Ex: 475nm, Em: 500–550nm) Green (Ex: 525nm, Em: 580–640nm) Red (Ex: 625nm, Em: 660–720nm) AFC (Ex: 405nm, Em: 495–505nm)
Detection Limit	2fmol fluorescein/200µl
Linear Dynamic Range	>6 logs (assay dependent)

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UV-Visible Absorbance Module Specifications

Detector	Top-reading, head-on photon counting photomultiplier tube (PMT)
Light Source	Xenon flash lamp
Spectra Range	200–600nm
Filter Wheel	GloMax [®] Discover: 9 absorbance filters included GloMax [®] Explorer: 5 absorbance filters
Wavelengths included	GloMax [®] Discover: 230, 260, 280, 320, 405, 450, 490, 560 and 600nm 10nm bandwidth GloMax [®] Explorer: 405, 450, 490, 560 and 600nm 10nm bandwidth
Detection Limit	0.1 O.D.
Dynamic Range	0–4.0 O.D.

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Product	Size	Cat.#
GloMax [®] Discover System <i>Configured with Luminescence, Fluorescence and UV-Visible Absorbance</i>	1 each	GM3000*
GloMax [®] Explorer System <i>Configured with Luminescence, Fluorescence and Absorbance</i>	1 each	GM3500*
GloMax [®] Explorer System <i>Configured with Luminescence and Fluorescence</i>	1 each	GM3510*

*Comes with U.S. Instrument and Tablet PC power cords.

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Product	Size	Cat.#
GloMax [®] Dual Injectors with Pump Station	1 each	GM3030
GloMax [®] Discover Luminescence Filter Paddle	1 each	GM3011
GloMax [®] Discover Fluorescence Filter Paddle	1 each	GM3012
GloMax [®] Discover Inner and Outer Injector Tubing, Natural	1 each	GM3014
GloMax [®] Discover Inner and Outer Injector Tubing, Black	1 each	GM3015
GloMax [®] Explorer Visible Absorbance Module Upgrade	1 each	GM3520
GloMax [®] Explorer GM3500 Upgrade to GloMax [®] Discover GM3000 <i>(only available for GM3500 model)</i>	1 each	GM3530
GloMax [®] Explorer GM3510 Upgrade to GloMax [®] Discover GM3000 <i>(only available for GM3510 model)</i>	1 each	GM3540

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

Caspase-Glo, CellTiter-Glo, Dual-Glo, GloMax, Kinase-Glo, Nano-Glo, ONE-Glo and P450-Glo are registered trademarks of Promega Corporation. ADP-Glo, CellTox, Chroma-Glo, DLR and NanoBRET are trademarks of Promega Corporation.

Windows is a registered trademark of Microsoft Corporation.

PROMEGA CORPORATION • 2800 WOODS HOLLOW ROAD • MADISON, WI 53711-5399 USA • TELEPHONE 608-274-4330

www.promega.com • © 2015 PROMEGA CORPORATION • ALL RIGHTS RESERVED • PRICES AND SPECIFICATIONS SUBJECT TO CHANGE WITHOUT PRIOR NOTICE • PUBLISHED IN USA 3/15 • 16342722 • PART BR261